

FOR EXPORT ONLY

ElaDerm

Microneedle RF



BISON
www.bisonmedical.com

The New Paradigm in RF Microneedling

**No pain microneedling RF with
Body Safe Cryo GAS**

ElaDerm Advantage

3 Key Benefits of Cryo Gas(BISON Patent)

1 No Pain

No more suffer from Pain in needle RF. ElaDerm works without pain with Body Safe Cryo Gas.

2 No Side Effect

ElaDerm causes no side effect such as burn, PIH, and scar due to the special Cryo Gas.

3 Enhanced Results

ElaDerm show enhanced and better results by using higher power without pain and side effect due to the special Cryo Gas.

☒ Vacuum Suction


Reducing Pain and increasing RF effect on thin part of skin such as forehead

☒ Smart Energy Control

You can use 1MHz, 2Mhz alone or combination of 1MHz and 2MHz to get better results.

☒ Pulse Mode

Minimizing damage of skin but delivering higher energy safely




New
Needle RF Handpiece



Why ElaDerm?

ElaDerm operates without Pain from Microneedle due to Special Cryo Gas which emits at the same time with microneedling. The microneedles are driven directly into the skin at a certain depth, and radiofrequency energy is then released inside of the skin. This gives a more uniform deep tissue warming that will then stimulate the collagen and new collagen fibers are produced. The needle tip contains very small gold needles acting as an internal conductor

How to do **PainFree** MNRF



Cooling gas release to targeted area together with microneedle

Indication

- Tightening Skin
- Lifting Skin
- Revitalizing Skin
- Improving Blemishes
- Reducing Wrinkle & Pore
- Reducing Pore
- Reducing Scars such as Acne scars
- Eliminating Osmidrosis



Specification

Microneedle RF System

Tip	Insulated and non-insulated tips
	 5x5  Bipolar RF needle
Dimension (WxDxH)	W 406 X D 578 X H 1644 (mm)
RF type	Bipolar needle RF (Radiofrequency)
RF Frequency	1MHz, 2MHz
RF Power	Up to 50Watt
Cryo cooling (duration)	Pre & Post spray (0~99ms)
Vacuum suction	Available while procedure

This information can be updated according to the product upgrade.

BISON

R&D 213 Office 703 Factory 903,904 (Gasan-dong, Ace High-end Tower 6), 234, Beotkkot-ro, Geumcheon-gu, Seoul, Korea (zipcode 08513) Tel)+82-2-865-7121 | sales@bisonmedical.com | www.bisonmedical.com

Copyright© BISON CO.,LTD. All rights reserved.

BSBEP23A